



TS-DDV SERIES

TABLE TOP STEAM
STERILIZER

B Series Table Top Steam Sterilizer(Automatic door)



Features

- ☆ Automatic door switch
- ☆ Automatic water
- ☆ Quality stainless steel inner tank,sealed cover
- ☆ Microcomputer control
- ☆ Safety interlocking device
- ☆ With drying function
- ☆ With BD test,HELIX test,Vacuum test.
- ☆ LED display and inductive operation
- ☆ Equipped with standard test interface
- ☆ European CLASS B standard.
- ☆ With automatic protection function: over temperature protection,over pressure protection, low water level protection, anti dry burning, safety door lock.
- ☆ Three times per-vacuum
- ☆ Buzzer reminder after sterilization,automatic stop.
- ☆ Automatic discharge of cold air,automatic steam exhaust after sterilization.
- ☆ Built-in steam generator to provide saturated steam quickly
- ☆ With built-in printer,USB function and U disk

Mode			
l Technical parameter	TS-DDV 18L	TS-DDV 23L	TS-DDV 45L
Volume	18L	23L	45L
Power	2kw	2kw	2.5kw
Voltage	110V/220V,60Hz/50Hz		
Design pressure	0.25MPa		
Design temperature	139℃		
Rated working pressure	0.22MPa		
Rated working temperature	134℃		
Chamber volume	0250×370	0250×470	0316×615
Outer dimension	620×505×450	505×710×450	657×870×575
Packing size(L×W×H)	670×560×510	760×560×510	980×750×630
Gross/net weight	62/53kg	65/58kg	142/105kg