

VA-SL (Digital Displayed Automatic Type)

Features

- ☆ Hand wheel translation and quick opening structure
- ☆ High quality stainless steel material
- ☆ Microcomputer automatic control, arbitrary setting sterilization parameters.
- ☆ Safety interlocking device
- ☆ Digital display and Touch buttons
- ☆ Standard test interface
- ☆ Self-expanding seal
- ☆ With automatic protection function: over temperature protection; over protection ; low water level protection; anti dry burning.
- ☆ With stainless steel bucket or basket
- ☆ Buzzer reminder after sterilization, automatic stop.
- ☆ Automatic discharge of cold air, automatic steam exhaust after sterilization.
- *Optional external printer



Model	VA-SL35	VA-SL50	VA-SL75	VA-SL100
Technical parameter				
Volume	35L	50L	75L	100L
Power	3.5kW	3.5kW	4.5kW	4.5kW
Voltage	AC 220V、50HZ			
Rated working pressure	0.25Mpa			
Rated working temperature	134℃			
Sterilization temperature	116℃~134℃			
Sterilization time setting range	4~120 min			
Chamber volume	Ø350×400	Ø350×525	Ø400×625	Ø450×650
Bucket dimension	Ø330×320	Ø330×480	Ø380×540	Ø430×550
Basket dimension (optional)	Ø320×350	Ø320×240×2	Ø360×280×2	Ø410×300×2
Outer dimension	480×480×1020	480×480×1050	525×525×1200	555×555×1250
Packing size (L× W× H)mm	580×630×1130	580×630×1230	630×680×1290	690×650×1410
Gross/net weight	79/58KG	106/80KG	120/90KG	133/106KG