

VERTICAL PULSE VACUUM STEAM STERILIZER

MODEL VA-SM



PRODUCT CATALOG

FEATURES

- Hand wheel translation and quick opening structure;
- **THREE TIME PRE-VACUUM;**
- High quality stainless steel material
- Microcomputer automatic control, arbitrary setting sterilization parameters;
- Safety interlocking device;
- LED touch screen display;
- Standard test interface;
- Self-expanding seal;
- **With drying function;**
- With stainless steel bucket or basket;
- With automatic protection function;
 - Over temperature protection;
 - Over pressure self discharge protection;
 - Low water level protection;
 - Anti drying burning
- Buzzer reminder after sterilization, automatic stop;
- **Automatic water;**
- Automatic discharge of cold air, automatic steam exhaust after sterilization, internal steam cycle;
- **External steam generator;**
- **Internal printer;**

技术参数 Technical parameter	VA-SM35	VA-SM50	VA-SM75	VA-SM100
容积 Volume	35L	50L	75L	100L
功率 Power	3.5KW	3.5KW	4.5KW	4.5KW
电压 Voltage	AC 220V、50HZ			
设计压力 Design pressure	0.25MPa			
设计温度 Design temperature	139℃			
额定工作压力 Rated working pressure	0.22MPa			
额定工作温度 Rated working temperature	134℃			
灭菌温度选择范围 Sterization temperature setting range	116℃-134℃			
灭菌时间选择范围 Sterilization time setting range	4~120 min			
干燥时间选择范围 Drying time setting range	0~240 min			
灭菌室容积 Chamber volume	Ø 350x400	Ø 400x505	Ø 400x655	Ø 450x650
内筒尺寸 Bucket dimension	Ø 330x320	Ø 380x430	Ø 380x560	Ø 420x540
(选配)网篮尺寸 Basket dimension(mm)	Ø 320x350	Ø 360x280x1	Ø 360x280x2	Ø 410x300x2
外形尺寸 Outer dimension(mm)	698x498x940	630x500x1180	725x525x1170	760x560x1260
包装尺寸 Packing size(LxWxH (mm))	790x610x1130	800x600x1230	840x610x1310	870x650x1380
毛/净 重 Gross/net weight	130/115 kg	138/118kg	150/124kg	159/132kg

