

# Ethylene Oxide Sterilizer

## SD-XH(120L, 220L)

---



### Features:

1. Sterilization chamber:5052 Aluminium, has better thermal conductivity than stainless steel
2. Vacuum system:High-performance, low-noise, oil-free vacuum pump
3. Built-in printer:Thermal printer
4. Humidity system:Use saturated steam to humidify the sterilization chamber, and the air humidity is more uniform
5. Timing system:EO gas diffusion time and residual EO aeration timing settable
6. Pre-set: Standard Sterilization Mode and Residual Washing

### Specifications:

- Sterilizing temperature:36-55 C adjustable
- Sterilizing Time:8-10 Hours , depends on materials
- Sterilization Method:Low temperature , gas sterilization .
- Display:LCD screen , Soft touch operate key
- Safety Features:Temperature and Pressure protection
- Tray:2 stainless net tray
- Voltage:220V, 50Hz
- Data Review:Thermal mini-printer and USB port (download)
- EO gas:100% Pure EO gas cartridge