

TS-A SERIES

TABLE TOP STEAM STERILIZER

A Series Table Top Steam Sterilizer



TS-A18/TS-A23

Features

- ☆ Quality stainless steel inner tank,sealed cover;
- ☆ The microcomputer controls the automatic cycle,121℃ and 134℃ can be set;
- ☆ Touch buttons;
- ☆ Equipped with standard test interface;
- ☆ European CLASSN standard;
- ☆ With automatic protection function:over temperature protection,over pressure protection,low water level protection,anti dry burning;
- ☆ Buzzer reminder after sterilization,automatic stop;
- ☆ Automatic discharge of cold air,automatic steam exhaust after sterilization;
- ☆ Built-in steam generator to provide saturated steam quickly;

Technical parameter	型号 Model	TS-A18	TS-A23
Volume		18L	23L
Power		1.5kw	2kw
Voltage		110V 220V60Hz / 50Hz	
Design pressure		0.25MPa	
Design temperature		139℃	
Rated working pressure		0.22MPa	
Rated working temperature		134℃	
Chamber volume		∅ 250×370	∅ 250×452
Outer dimension		500×400×350	600×400×350
Packing size (L×W×H)		540×450×390	640×450×390
Gross/net weight		26/24kg	29/27kg