

Table Top Pressure Steam Autoclave

TS-B



Features:

1. Sterilizing course: Automatic sterilization controlled by computer, easy to operate.
2. Maximum temperature: up to 134° C, suit for 4-6 minutes rapidly sterilizing.
3. Safe protection of water lacking.
4. With over-temperature, over-pressure auto-protect device:
5. After sterilization it will automatically cut the power and then alarm.
6. Door safety lock.
7. Fully stainless steel structure.

Specifications:

- Capacity: 16L, 24L
- Working pressure: 0.22Mpa
- Working temperature: 134°C
- Adjustment of temperature: 105-134°C
- Timer: 0-60 min
- Power: 1.5KW

