

Table Top Pressure Steam Autoclave

TS-CD(Digital display)



Features:

1. For 4~6 minutes rapidly sterilizing.
2. Digital display of working status, touch key
3. With 3 fixed sterilizing programs and one adjustable program for user;
4. The whole process of water adding, temperature rising, sterilizing, drying steam discharging controlled automatically.
5. Steam-water inner circulation system: no steam discharge, and the environment for sterilizing will be clean and dry.
6. Exhaust the cool air automatically.
7. Safe protection of water lacking.
8. Door safety lock system.
9. With four stainless steel sterilizing plates.
10. The chamber of the sterilizer is made of stainless steel.
11. Automatically shut off with beep reminding after sterilization finished.
12. With drying function.

Specifications:

- Capacity: 20L, 24L, 35L, 50L
- Working pressure: 0.22Mpa
- Working temperature: 134°C
- Adjustment of temperature: 105-134°C
- Timer: 0-99 min
- Power: 1.5-2.7KW

