

Table Top Pressure Steam Autoclave

TS-DDV(Class B)

Features:

1. European Standard Class B with three times vacuum and drying, the remaining temperature of sterilized instrument is less than 0.2%.
2. LCD display working status, touch keys.
3. The vacuum data can reach -0.8bar, it fits various instrument, including package, unpackaged, solid, hollow, multi-apertures and inset pipes. It makes sure the instrument can be sterilized completely.
4. Fully computer controlled LCD display bright show. The interface is modularized panel for easy operation.
5. It is installed with BOWIE & DICK which measures the penetration of water steam.
6. It has vacuum test which can make sure the ability of vacuum testing.
7. Built-in fast and independent steam generator.
8. Built-in mini printer which can record the sterilizer information.
9. Safety valve, Safety lock door system, Pressure and temperature overload protection, Low water protection, Over current protection, Alarm system.
10. Can inactivate HIV, HBV, BSE, and bacillus efficiently.

Specifications:

- Capacity: 18L, 24L
- Working pressure: 0.22Mpa
- Working temperature: 134°C
- Adjustment of temperature: 105-134°C
- Timer: 0-99 min
- Power: 2KW

