

Vertical Pressure Steam Autoclave

VA-FJ(Indicator Light Display)



Features:

1. Indicator light indicates working state
2. Over-temperature & over-pressure auto-protection device
3. Double scale indication pressure gauge
4. Safe protection of water lacking
5. With the below exhaust structure, which is convenient to exhaust the cold air
6. Automatically shut off with beep reminding after sterilization
7. Fully stainless steel structure
8. Easy to operate, safe and reliable
9. With two stainless steel sterilizing baskets

Specifications:

- Capacity: 35L, 50L, 75L, 100L, 150L
- Max. Working pressure: 0.23Mpa
- Working pressure: 0.22Mpa
- Working temperature: 134°C
- Adjustment of temperature: 105~134°C
- Timer: 0-60min
- Power: 2.5-6KW