



Hefei Sada Medical Equipment Co., Ltd.
Address: No. 12629, Jingdian Plaza, Fanhua Avenue,
Economic Development Zone, Hefei, Anhui, China
Tel/Fax: +86-551-63879990
Website: www.sadamedical.com

VA-FB (Digital Displayed Automatic Type)



Features

- ☆ Bolt structure
- ☆ High quality stainless steel material
- ☆ Microcomputer automatic control, arbitrary set of sterilization parameters
- ☆ Safety interlocking device
- ☆ Digital display and touch buttons
- ☆ Standard test interface
- ☆ Self-expanding seal
- ☆ The container with automatic protection function: over temperature protection; over pressure protection self discharge; low water level protection, anti dry burning
- ☆ With stainless steel bucket or basket
- ☆ Buzzer reminder after sterilization, automatic stop.
- ☆ Automatic cooling air, sterilization end automatic exhaust steam

*Optional external printer

Model	VA-FB35	VA-FB50	VA-FB75	VA-FB100	VA-FB120	VA-FB150
Technical						

parameter						
Volume	35L	50L	75L	100L	120L	150L
Power	3.5kW	3.5kW	4.5kW	4.5kW	4.5kW	6.0KW
Voltage	AC 220V、50HZ					
Rated working pressure	0.22Mpa					
Rated working temperature	134℃					
Sterilization temperature	116℃~134℃					
Sterilization time setting range	4~120 min					
Chamber volume	Ø300×500	Ø350×525	Ø400×625	Ø450×650	Ø450×760	Ø500×780
Bucket dimension	Ø280×400	Ø330×480	Ø380×560	Ø420×540	Ø420×660	Ø470×720
Basket dimension (optional)	Not optional	Ø320×240×2	Ø360×280×2	Ø410×300×2	Ø410×350×2	Ø460×350×2
Outer dimension	400×400×860	480×480×960	525×567×1020	555×555×1080	555×599×1200	605×649×1210
Packing size (L× W× H)mm	550×500×1030	580×630×1100	630×680×1190	690×650×1230	690×650×1330	750×700×1370
Gross/net weight	58/41KG	75/53KG	82/61KG	103/78KG	108/82KG	116/90KG