

Vertical Pressure Steam Autoclave

MODEL: VA-FJ



Technology Advantage:

1. Fully stainless steel 304 structure.
2. Microprocessor control, indicator light indicates working status: pressure, time, temperature.
3. Auto discharge the cool air, and steam discharging automatically after sterilizing.
4. Double scale indication pressure gauge.
5. Over temperature and over pressure auto-protection.
6. Safe protection of water lacking.
7. Self-inflating type seal.
8. Automatically shut off with beep reminding after sterilization.
9. With two stainless steel sterilizing baskets.

Safety Device:

1. Safe protection of water lacking
2. The program will be stopped if water level lower than the low limit.
3. The buzzer alarms when a cycle of operation is completed.
4. Over temperature&over pressure auto-protection.
5. Over current protector.
6. Safety valve of auto release against over-pressure.

Accessory:

1. Stainless steel basket
2. Fuse
3. Hose
4. Seal ring
5. Wrench

Main Parameters:

Model	VA-FJ35	VA-FJ50	VA-FJ75	VA-FJ100	VA-FJ120	VA-FJ150
Volume	35 Liters	50 Liters	75 Liters	100 Liters	120 Liters	150 Liters
Sterilizing Chamber Volume	Ø318*450mm	Ø340*550mm	Ø400*600mm	Ø440*650mm	Ø480*660mm	Ø510*740mm
Working Pressure	0.22 Mpa					
Max Working Pressure	0.23 Mpa					
Working Temperature	134°C					
Adjustment of Temperature	105°C-134°C					
Timer	0-60 minutes					
Heat Average	≤± 1°C					
Power	2.5 KW AC220V 50Hz	3 KW AC220V 50Hz	4.5 KW AC220V 50Hz		6 KW AC220V 50Hz	
Overall Dimension(mm)	480*460*850	520*520*980	560*560*980	590*590*1080	600*640*1140	670*690*1130
Transport Dimension(mm)	570*550*970	590*590*1110	650*630*1150	680*650*1220	730*730*1270	760*760*1270
Gross Weight/Net Weight	56kg/42kg	68kg/50kg	90kg/70kg	105kg/85kg	125kg/100kg	135kg/110kg